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取样：**Isolation Valve → Exchange → Start → Vent**

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确认 Beam：**Isolation Valve → Faraday Cup → Measure → Reference → Focus**

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### 1. Exposure preparation

Confirm the vacuum of SIP and Gun Monitor parameters. Record the vacuum of the instrument and register. Start SEM and WecaS software.

**Unload the Cassette: Isolation Valve→ Exchange→ Start→ Vent**

Close Isolation Valve and move the stage to the Exchange position by clicking Exchange;

Check the Isolation Valve sensor is lighting to confirm the closed state;

Confirm that the START button is blue, UNLOAD and STAGE are lit green. Press the blue START button; LOAD and EXCHG are lit. Press VENT button. The LED of VENT button turns off when the atmosphere pressure is reached;

Open the top cover and remove the cassette;

Push the EVAC button and pull the vacuum of the exchange room to finish.

**Load the Cassette: Mounting the sample→ VENT → Set the cassette→ START**

Mounting a sample onto the Cassette. Be careful not to scratch the Cassette;

Clean the sample by N<sub>2</sub> gun;

Press VENT button. The LED of VENT button turns off when the atmosphere pressure is reached;

Open the top cover;

Set the cassette. Make sure that the cassette is not shifted by moving the cassette to the right or left;

Confirm that the START button is blue, LOAD and EXCHG are lit green. Press the blue START button; UNLOAD and STAGE are lit. The color of START button is blue, loading is completed.

**Confirm Beam Condition: Isolation Valve→ Faraday Cup→ Measure→ Reference→ Focus**

Click Faraday Cup button to move the stage to the Faraday Cup position;

Open the Isolation Valve if the pressure of the main Chamber is better than 1E-4 Pa;

Click Measure button to measure beam current and adjust beam current by click the arrow;

Move the stage to the Reference position;

Adjusting the Focus and Stigmator;

### 2. Parameter Settings

The SCON file is necessary for the exposure.

**SCON File preparation: Start Beamer software→ Import the gds file→ Edit the pattern→ Export SCON file→ Confirm the Dot, Field size and Pitch values**

### 3. Exposure

The exposure procedure can be adjusted as required.

**The Exposure Procedure: Start WecaS software→ Edit Schedule→ Position→ Calculate Dose Time→ Set parameters→ Select Sample Holder→ Move pattern→ Estimated Time→ Expose→ Close Isolation Valve.**